

# GemmaCert HPLC Protocol: Analytical Method for Determination of Cannabinoids



**1. Equipment, Solvents, Reagents and Devices:**

**1.1. Solvents, Reagents, Analytical standards and disposable Equipment:**

Methanol HPLC grade  
Acetonitrile HPLC grade  
Formic Acid 99%  
Water HPLC grade  
Ammonium formate  
Sigma (Cerilliant) /Restek analytical standard  
50 mL centrifuge tubes  
15 mL centrifuge tubes  
PTFE Syringe Filter 0.22  $\mu\text{m}$   
0.45  $\mu\text{m}$  *membrane filter* .  
Clear and amber glass HPLC vials

**1.2. Devices:**

Sonicator  
Shaker  
MERCK 'HITACHI' Lachrome HPLC system (pumps, injector, UV detector)/  
MERCK 'HITACHI' Lachrome Elite HPLC system (pumps, autosampler, UV detector)

**2. Chromatographic Conditions:**

Column: Raptor ARC-18 2.7  $\mu\text{m}$  150x4.6mm

Guard- Raptor ARC-18 2.7  $\mu\text{m}$  5x4.6mm

Detector: UV at 228nm

Column temp: 30°C

Autosampler temp: 10°C (N.A for Hitachi)

Flow rate: 1.5 ml/min

Flow mode: Isocratic A: B (25:75)

Injection volume: 5.0  $\mu\text{L}$

Run time: 10.0 min

Mobile phase-Preparation for 1L solution:

Mobile Phase A- 10mM Ammonium Formate and 0.1% Formic Acid Buffer:

Weigh about 0.6 gr of Ammonium Formate, dissolve with 1000ml water and add 1.0 ml of Formic Acid 99%-filter mobile phase A through 0.45  $\mu\text{m}$  *membrane filter* .

Mobile Phase B-0.1% Formic Acid in Acetonitrile:

Add 1.0 ml of Formic Acid 99% into 1000ml of Acetonitrile.

Mobile phases are stable if system suitability can be achieved and not more than 2 weeks at RT.

Diluent: Methanol

Injection wash: Methanol

Pressure: 0-300bar maximum.

**Table 2.1:**

Time (min)	%A	% B (ACN)	% C (10% ACN in Water)	% D (ACN)
0	25	75	0	0
10	25	75	0	0

**3. Standard preparation:**

Stock sol 1:

Accurately transfer 50 $\mu$ L of each standard (conc.: 1 mg/mL =1000ppm) into the same HPLC vial add 100 $\mu$ L diluent, Mix with vortex. This is stock solution 1 (conc. 100ppm). Dilute next standards with  $V_{\text{Stock solution}}$  and  $V_{\text{diluent}}$  according to following table, after each dilution mix vial with vortex:

**Table 4.1:**

Concentration [ppm]	$V_{\text{Stock solution (100ppm)}}$	$V_{\text{diluent}}$	Total volume:
100ppm stock sol 1			500 $\mu$ L
50ppm	100 $\mu$ L	100 $\mu$ L	200 $\mu$ L
25ppm	50 $\mu$ L	150 $\mu$ L	200 $\mu$ L
10ppm-stock solution no.2	150 $\mu$ L	600 + 750 $\mu$ L	1500 $\mu$ L
Use stock solution no.2 for next dilutions:			
5ppm	100 $\mu$ L	100 $\mu$ L	200 $\mu$ L
2.5ppm	50 $\mu$ L	150 $\mu$ L	200 $\mu$ L
1ppm	50 $\mu$ L	450 $\mu$ L	500 $\mu$ L

\*Any other combination of dilutions may be used if final conc. Is maintained.

**4. System suitability:**

- a. New 7<sup>1</sup> points calibration curve must be injected as system suitability before any new project unless QC solution (10 ppm) is injected as system suitability and Total %RSD when taking into consideration last QC solution from previous calibration curve and all following QC solution injections during last run is under 2.5% for all the 8 cannabinoids<sup>2</sup>.
  - 1- One point from the calibration curve can be removed if necessary (excluding 100 ppm) to achieve linearity.
  - 2- System suitability must be checked prior to samples injection, when %RSD is not met, a new calibration curve must be injected.
- b. Calibration's curve 10ppm solution will be injected in triplicate to demonstrate the system's reproducibility and stability. %RSD for all 8 cannabinoids must be NMT 2.5%.
- c. Retention time and order of elution of main peaks is according to next table:

**Table 4.1:**

No.	Cannabinoids Name	Approx. Retention time (min)
1	CBDA	2.7
2	CBGA	2.9
3	CBG	3.0
4	CBD	3.1
5	THCV	3.4
6	CBN	4.6
7	$\Delta^9$ -THC	5.8
8	THCA	7.6

- d. Tailing for all the peaks should be NMT 2.0 and NLT 0.8. Tailing may be examined for 10ppm and higher standard solutions.
- e.  $R^2$  for all cannabinoids must be NLT 0.999

**5. Sample preparation:**

5.1. Extraction: Grind the flower using suitable grinder with at least three rotations and transfer about 0.2 gr of grinded flower into a 50-mL centrifuge tube (For big flowers weigh in duplicates or triplicates). Dissolve the sample by Accurately adding 10.0 mL of diluent into the centrifuge tube. Sonicate for 15 min in iced water and then shake for another 15 min. Filter about 1.5mL of the sample solution into an amber HPLC vial through 0.22 PTFE membrane filter. This is stock sample solution. 0.02 gr/ml

5.2. dilutions:

5.2.1. Dilution factor 100: transfer 100  $\mu$ L of stock sample solution and 9.9 mL of diluent into a 15-mL centrifuge tube, mix with vortex and pipet a portion of the solution into a HPLC vial. 0.0002 gr/ml

5.2.2. Dilution factor 8: transfer into a HPLC Vial 130 $\mu$ L of stock sample solution with 910 $\mu$ L diluent and mix with vortex. 0.0025 gr/ml

5.3. Oil Extraction:

For full oil extraction process follow OIL00001 work instructions.

For analyses of oil extract, extract as described in section 6.1 then, dilute the extracts by factor 100 (as in section 6.2.1)-this is intermediate solution. Further dilute the samples by factor 5:

Accurately transfer 200  $\mu$ L of the intermediate solution + 800  $\mu$ L MeOH to an HPLC vial, Mix. 0.00004 gr/ml

**6. Sample for LOD test preparation:**

6.1. For each cannabis species prepare sample for LOD in triplicate. Weigh an empty weighing vessel, zero the balance and weigh about 2.0 gr of grind cannabis collected from each species flowers. Place the weighed vessels with samples, inside an oven heated to 70°C for 72 hours. (According to LOD00001 procedure)

6.2. Weigh the weighing vessels with dried sample and calculate %LOD:

$$\%LOD = \frac{W1 - W2}{Sample\ weight} * 100\%$$

W1 – Vessel weight before drying

W2 – Vessel weight after drying

6.3. %LOD can be also tested by using sartorius Moisture analyzer according to SOP00004 and LOD00001 procedures.

**7. HPLC PROCEDURE:**

- 7.1. Inject diluent at least twice to identify system peaks and any interferes with the method peaks of interest (appendix 1).
- 7.2. Inject calibration curve or QC solution and check system suitability according to section 5 (appendix 2).
- 7.3. Inject samples, if both concentrations are injected together the lower concentration samples should be injected first.
- 7.4. After every 20 samples inject diluent and then QC solution, Inject QC solution in the end of the sequence as well. % RSD for all QC injection must be NMT 2.5%<sup>1,2</sup>.
- 7.5. In the end of work use a proper program for washing the column to manufacture storage instructions.
  - 1- If % RSD was Higher than 2.5%, wash column with 100% ACN (1.5 ml/min). Then set again to A: B (25:75) mobile phase. Wash autosampler and inject MeOH in triplicate and one QC (10ppm solution). Check %RSD of this injection and 3 replicates from calibration curve. If %RSD is less than 2.5% work is valid or can be continued in case of injecting part 2.
  - 2- Removing one QC solution- if one of the QC solutions injected during the sequence (**excluding the one injected in the end**) is clearly exceptional (by peak area) than all the others, which demonstrate a clear trendline, this solution may be removed for %RSD calculations. This may occur due to pressure drop or some failure in injection process.

**8. Integration:**

- 8.1. Ignore any peak which is not injected in the calibration curve.
- 8.2. Integration may be inhibited before elution of CBDA and after the elution of THCA.
- 8.3. Minimum area should be no bigger than the area of the smallest peak in 1.0ppm standard solution injection.

**9. Calculations and reports:**

- 10.1 Calculate sample concentration (ng/μL) by validated calculating program or by HPLC software.

Calculate %W according to following equation:

$$\%W = \frac{\frac{Amount[ng]}{inj\ vol[\mu L]} \cdot Extr.\ vol[ml] * 10^3[ml \rightarrow \mu L] \cdot DF \cdot 100\%}{M * 10^9[gr \rightarrow ng] \cdot Purity}$$

*Extr. vol* – extraction vol [ml]

*DF* – dilution factor


*M* – dried sample weight [gr]  $M = Sample\ weight - \left(\frac{Sample\ weight \cdot \%LOD}{100}\right)$

*Purity*–standard potency/Assay

- 10.2. Report all relevant data to system suitability, chromatograms, sample concentration and %W for all samples.

Appendix 1:

Blank Chromatogram:



**Appendix 2:**
**100ppm Chromatogram:**

-for CBG Asymmetry see validation report -manual calculation.

<b>28 10STD 100ppm-3</b>					
Sample Name:	10STD 100ppm-3	Injection Volume:	5.0		
Vial Number:	RC7	Channel:	UV_VIS_1		
Sample Type:	unknown	Wavelength:	228		
Control Program:	Restek ARC-18 20_06_19	Bandwidth:	n.a.		
Quantif. Method:	Restek ARC-18_8STD	Dilution Factor:	1.0000		
Recording Time:	23/6/2019 15:32	Sample Weight:	1.0000		
Run Time (min):	10.01	Sample Amount:	1.0000		

Linearity\_Restek 8STD #28 [modified by Lenovo] UV\_VIS\_1  
WVL:228 nm

1 - CBDA - 2.660  
2 - CBGA - 2.853  
3 - CBG - 2.987  
4 - CBD - 3.133  
5 - THCV - 3.393  
6 - CBN - 4.607  
7 - THC-9 - 5.747  
8 - THCA - 7.520

No.	Ret.Time min	Peak Name	Height mAU	Area mAU*min	Rel.Area %
1	2.66	CBDA	422.653	23.507	16.22
2	2.85	CBGA	388.841	22.394	15.45
3	2.99	CBG	209.719	12.567	8.67
4	3.13	CBD	208.342	12.609	8.70
5	3.39	THCV	190.517	11.882	8.20
6	4.61	CBN	382.275	29.226	20.17
7	5.75	THC-9	123.810	11.297	7.80
8	7.52	THCA	173.945	21.420	14.78

No.	Ret.Time min	Area uAU*Sec	K'	Resolution(EP)	Asymmetry(EP)
1	2.66	1410401.336	2.06	2.19	1.34
2	2.85	1343662.618	2.28	1.46	1.34
3	2.99	754042.016	2.43	1.51	n.a.
4	3.13	756540.058	2.59	2.82	1.18
5	3.39	712901.962	2.90	11.33	1.28
6	4.61	1753574.540	4.30	8.79	1.20
7	5.75	677810.276	5.61	10.65	1.21
8	7.52	1285193.844	7.64	n.a.	1.44

Change control table:

Edition	Details of change	Effective date
1.0->2.0	Section 2, line 12 addition of 'ml' units	03.06.20
2.0 => 3.0	Section 1.2, changed Dionex to Lachrome Elite HPLC	19.10.20

Document approval table:

Position	Name	Signature
Lab Manager	Olga Pendo	
CEO	Guy Setton	